



SPECTRO-V11D Spectrophotometer

SPECTRO-V11D is the only model of manually setting wavelength among families, But precise design & high quality components ensures excellent performance. It is widely used in high schools and colleges for general analysis and experiments

- Optical System: Single Beam, Grating 1200 lines/mm
- Detector: Silicone Photo diode

Features: Large LCD Screen(128x64Dots) • The wavelength can be read out from the screen directly • Auto zero & blank, easy to use. Press one button for easy switching of Transmittance, Absorbance & Concentration modes • Parallel port, data can be exported to printer directly • Coefficient method to measure the

unknown samples. $C=K \times A \times B$. It can memorize the Last input of K & B • Large sample compartment, it can accommodate 5-100mm path length cuvettes with optional holders • Pre-aligned design ensures the user to change lamp conveniently • Optional software M.Wave Professional based on windows® can expand the applications to Standard Curve & Kinetics • High quality silicon photometric diode detector and 1200 lines/mm grating ensure the high accuracy and precision.

Model	SPECTRO-V11D	SPECTRO-UV11
Wavelength range	325-1000nm	200-1000nm
Spectral Bandwidth	4nm	
Wavelength Accuracy	±2nm	
Wavelength Repeatability	1nm	
Wavelength Setting	Manual	Auto
Photometric Accuracy	±0.5%T	
Photometric Repeatability	0.3%T	
Photometric Range	-0.3-A, 0-200%T	
Stray Light	0.3%T	
Stability	±0.004A/h @ 500nm	
Display	128 x 64 Dots LCD	
Photometric Mode	T, A, C, F	
Standard Cell Holder	4-position 10mm cell changer	
Sample Compartment	Standard 10mm path length cuvette	
Light Source	Tungsten Lamp	Tungsten&Deuterium lamp
Output	USB Port&Parallel Port (printer)	
Power Requirement	AC 220V/50Hz or AC110V/60Hz	
Dimensions (mm)	W480xD360xH160	
Weight	10kg	12kg


Wave Professional PC-Control Software

The windows® based PC application software takes the best features of the stand-alone version plus more powerful data processing and expanded data collecting and storage capability. The software is optional. Once it be used, the uv/v-l series become more excellent



Quantitative: Use up to 38 standards to establish standard curve. Four methods for fitting a curve: Linear fit, Linear through zero, Square fit.

Kinetics: The Kinetics mode may be used for time course scanning or reaction rate calculations. Abs. Vs. time graphs is displayed on the screen in real time. Wait time, test time and time intervals may be set.



Wavelength Scan:

Automatically record peaks & valleys. The quantity of the curves is unlimited that you store.

Post- run manipulation & processing includes: Re-scaling axes, curve Smoothing, combination, zooming, overlap... 1 st to 4th derivative.

Multi-wavelength Test: You can set at most 20 wavelengths to measure a sample.

DNA Test:

2 formulas are under your choice:

ONACon. = $62.9 \times A_{260} - 36.0 \times A_{280}$
 or = $49.1 \times A_{260} - 3.48 \times A_{230}$
 Protein Con. = $1552 \times A_{260} - 757.3 \times A_{280}$
 or = $183 \times A_{260} - 75.8 \times A_{230}$
 You can also enter other wavelengths and factors to calculate.



SPECTRO-UV/V12

SPECTRO-UV/V12, Spectrophotometer

The UV/V12 spectrophotometer have the most attractive cost performance. Using your standard sample solutions, you can get a standard curve on the large LCD screen by the local control software, also you can print the curve through the parallel port. They are widely used in colleges and enterprises for general quantitative analysis and experiments.

Optical System: Single Beam, Grating 1200 lines/mm • Detector: Silicone Photo diode • SPECTRO-V12 Range: 325-1000nm, SPECTRO-UV12 Range: 200-1000nm.

Functions:

1. Basic Mode Absorbance, Transmittance or Concentration measurements.

2. Quantitative I) **Standard Curve** At most 9 standard samples can be used to establish a standard curve, the curve and the curve equation will display on the screen simultaneously. You can measure your unknown concentration solutions by the curve. II) **Coefficient Method** If you have known the coefficient k&b of the formula $C=kA+b$, you can input the value directly by the button, Then to test your unknown solutions.

Model	SPECTRO-V12	SPECTRO-UV12
Spectral Bandwidth	4nm	
Wavelength Accuracy	±2nm	
Wavelength Repeatability	0.8nm	
Photometric Accuracy	±0.5%T	
Photometric Repeatability	0.3%T	
Photometric Range	-0.3-3.0A, 0-200%T	
Stray Light	0.3%T	
Stability	±0.002A/h @ 500nm	
Sample Compartment	Standard 10mm pathlength cuvette	
Light Source	Tungsten Lamp	Tungsten&Deuterium lamp
Output	USB Port&Parallel Port (printer)	
Power Requirement	AC 220V/50Hz or AC110V/60Hz	
Dimensions (mm)	W470xD370xH180	
Weight	12kg	12kg



SPECTRO-UV3 Series

SPECTRO-UV3, Scanning Spectrophotometer

SPECTRO-UV3 Series is an advanced single beam design consisting of 10 models. They differ In bandwidth & wavelength accuracy, but provide excellent performance for measurements in the range of 190nm to 1100nm. They are suitable for clinical lab applications, pharmaceutical, and bio-chemical. as well as routine applications such as Quantitative analysis, Kinetics, Wavelength Scan, Multi-Wavelength, and DNA/Protein analysis. UV-Vis Analyst application software based Microsoft Windows makes these instruments versatile. All instruments provide excellent performance for measurements. They are divided into in two types: PC models and stand-alone models -In Stand-alone models, all

software methods are included as built-in standard; this eliminates the need of software. Online software update via Internet Data can be downloaded . The PC models come standard with Windows® based application software UV-Vis Analyst.

Features: Fixed or variable slits (bandwidths) Sealed, solvent-resistant tactile keypad with alpha-numeric entry for file names and units. Pre-aligned deuterium lamp for easy lamp replacement The status of the lamps may be monitored. Powerful built-in program or PC Windows™ based software UV/VIS. Analyst including sophisticated utility programs. Data Download-to-PC software for stand-alone models (optional). Real-time clock for date and time stamping of results.

Model	SPECTRO-UV30/UV30PC	SPECTRO-UV31/UV31PC	SPECTRO-UV32/UV32PC	SPECTRO-UV32S/UV32PCS	SPECTRO-UV33/UV33PC
Wavelength range	190-1100nm				
Spectral bandwidth	4nm	2nm	1.8nm	0.5/1/2/4nm	1nm
Optical system	Single beam, grating 1200 lines/mm				
Wavelength accuracy	±0.5nm		±0.3nm		
Wavelength repeatability	0.3nm		0.2nm		
Scan speed	Hi, Med, Low, Max.3000nm/min				
Photometric accuracy	≤±0.5%T or ±0.003A@1A				
Photometric range	-0.3 -3A, 0-200%T				
Stray light	≤0.05%T@220nm, 360nm				
Stability	±0.002A/h@500nm				
Display	5 inches LCD (320x240 dots)				
Baseline flatness	±0.002A(200-1000nm)				
Standard cell holder	Standard 10mm pathlength cuvette				
Light source	Halogen&Deuterium lamp (pre-aligned)				
Output	USB Port¶llel port (printer)				
Power requirement	AC 110/220V 50/60Hz				
Dimensions (mm)	W480xD360xH160			W600xD450xH200	
Weight	14kg			20kg	



SPECTRO-16/18 Series, Spectrophotometer

The SPECTRO-16/18 Series are developed for precise test, its stray light is only 0.05%T. They are flexible, easy to use and maximize value. The value is evident from the performance provided every day. The local stand-alone software provide functions of Basic Mode, Quantitative Test, Kinetics and System Utilities.

• Optical System: Single Beam, Grating 1200 lines/mm • SPECTRO-V16/V18 Range: 320-1100nm, SPECTRO-UV16/UV18 Range: 190-1100nm.

Features:

Large LCD screen (128x64Dots). Can display total 200 groups of data, 5 groups per screen. Can display standard curve and kinetics curve • System can also

save the test results, total 200 groups of data and 100 standard curves can be saved in the RAM. Be Convenient for check and reload • Data can be restored after a sudden power cut • Auto setting wavelength • Tungsten lamp & deuterium lamp can be turned on/off individually to extend lifetime • The optional application software M.Wave Professional provides complete control of the spectrophotometer from a Computer through the Built-in USB port. It can expand to the following functions: Quantitative, Kinetics, Wavelength Scan, Multi-wavelength & DNA/Protein Test • Pre-aligned design makes it convenient to Change lamps • Large sample compartment, It can accommodate 5-100mm path length cuvettes with optional holders. A variety of optional holders.

Model	SPECTRO-V16	SPECTRO-V18	SPECTRO-UV16	SPECTRO-UV18
Spectral Bandwidth	4nm	2nm	4nm	2nm
Wavelength Accuracy	±0.5nm			
Wavelength Repeatability	0.3nm			
Photometric Accuracy	±0.3%T			
Photometric Repeatability	0.2%T			
Photometric Range	-0.3-3A, 0-200%T			
Stray Light	0.05%T @ 360nm		0.05%T @ 220nm, 340nm	
Stability	±0.002A/h @ 500nm			
Display	Graphic LCD (128 x 64 Dots)			
Keyboard	22 Membrane keypad			
Sample Compartment	Standard 10mm pathlength cuvette			
Light Source	Tungsten Lamp		Tungsten&Deuterium lamp	
Output	USB Port&Parallel Port (printer)			
Power Requirement	AC 220V/50Hz or AC110V/60Hz			
Dimensions (mm)	W470xD370xH180			
Weight	14kg		14kg	

Accessories - Holder:

P.N 900410  4-Cell Holder for 10mm Squ.Cuvette	P.N 900420  4-Cell Holder for up to 50mm Squ.Cuvette	P.N 900430  4-Cell Holder for up to 100mm Squ.Cuvette	P.N 900540  Cylindrical Cell Holder	P.N 900610  Water-Jacketed Single Cell Holder
P.N 900210  Micro Cell Holder	P.N 900530  Test Tube Holder	P.N 900310  8-Position Auto Cell Changer	P.N 900550  Solid Sample Holder (Single Cell)	P.N 900620  10mm Water-Jacket 4-Cell Holder

Accessories - Others:

P.N 900910  Thermal Printer	P.N 900920  Stylus Printer	1. P.N 900120 2. P.N 900110 3. P.N 900130  1.Peltier Unit 2.Ambier Sipper Unit 3.Peltier/Sipper System
---	--	--

SPECTRO-UV6 Series, Double Beam Spectrophotometer

UV-6 Series are advanced double beam spectrophotometer consisting of four models:

Stand-alone models: SPECTRO-UV61

with 1.8nm fixed bandwidth & SPECTRO-UV63

with 1.0nm fixed bandwidth;

PC models: UV61PC with 1.8nm fixed

bandwidth & UV63PC with 1.0nm fixed

bandwidth. Other specifications of the four

models are almost the same except bandwidth. The two detectors are measuring sample and reference respectively and simultaneously for optimizing measurement accuracy. They provide excellent performance for measurements in the range of 190nm to 1100nm. They are suitable for pharmaceutical, biochemical and clinical lab applications as well as routine applications such as quantitative analyses, kinetics, spectrum scanning, multiple components and DNA/Protein, PC Windows application software make these instruments versatile. All instruments provide excellent performance for measurements.


- Optical System: Single Beam, Grating 1200 lines/mm

- To Stand-alone models, All software methods are

- included as built-in standard, thus eliminating the need

- for software options • Online software upgrade via

- internet helps to keep your software up-to-date •



- Data Download-to-PC software expands the data storage to unlimited • The stand-alone models with 5 inches screen and the PC models with 3.8 inches screen • Model SPECTRO-UV61/UV61PC Range: 190-1100nm • Model SPECTRO-UV63/UV63PC Range: 190-1100nm.

Model	SPECTRO-UV61	SPECTRO-UV61PC	SPECTRO-UV63	SPECTRO-UV63PC
Spectral Bandwidth	1.8nm	1.8nm	1.0nm	1.0nm
Wavelength Accuracy	±0.3nm			
Wavelength Repeatability	0.2nm			
Scanning Speed	Hi, Med, Low, Max, 3000nm/m in			
Photometric Accuracy	±0.3%T			
Photometric Repeatability	0.15%T			
Photometric Range	-0.3-3A, 0-200%T, 0-9999Conc.			
Stray Light	0.05%T			
Stability	±0.001A/h @ 500nm			
Display	5 inches LCD (320x240)	3.8 inches LCD (320x240)	5 inches LCD (320x240)	3.8 inches LCD (320x240)
Baseline Flatness	±0.001A(200-1000nm)			
Sample Compartment	Accommodates 100mm pathlength cuvette with optional holder			
Light Source	Halogen&Deuterium lamp(Pre-aligned)			
Output	USB Port&Parallel Port (printer)			
Power Requirement	AC 220V/50Hz or AC110V/60Hz			
Dimensions (mm)	W600xD450xH200			
Weight	22kg		22kg	


The PC models come standard with windows® based application software UV/VisAnalyst.

Accessories - Cells:

Square Cuvette		
Name&Specifications	P.N	Remark
Square Cuvettes.Glass 10mm	916101	/4pcs
Square Cuvettes.Glass 20mm	916102	/4pcs
Square Cuvettes.Glass 30mm	916103	/4pcs
Square Cuvettes.Glass 50mm	916104	/4pcs
Square Cuvettes.Glass 100mm	916105	/4pcs
Square Cuvettes.Quartz 10mm	916111	/2pcs
Square Cuvettes.Quartz 20mm	916112	/2pcs
Square Cuvettes.Quartz 30mm	916113	/2pcs
Square Cuvettes.Quartz 50mm	916114	/2pcs
Square Cuvettes.Quartz 100mm	916115	/2pcs

Micro Cell		
Name&Specifications	P.N	Remark
100UL Micro Cell	916121	/1pcs
Square Cuvettes.Glass 20mm	916122	/1pcs
Square Cuvettes.Glass 30mm	916123	/1pcs

Flow Cell		
Name&Specifications	P.N	Remark
5mm SEIF Masking Cont. Flowthrough G.Cell	916131	/1pcs
10mm SEIF Masking Cont. Flowthrough G.Cell	916132	/1pcs
20mm SEIF Masking Cont. Flowthrough G.Cell	916133	/1pcs
30mm SEIF Masking Cont. Flowthrough G.Cell	916134	/1pcs
5mm SEIF Masking Cont. Flowthrough Q.Cell	916141	/1pcs
10mm SEIF Masking Cont. Flowthrough Q.Cell	916142	/1pcs
20mm SEIF Masking Cont. Flowthrough Q.Cell	916143	/1pcs
30mm SEIF Masking Cont. Flowthrough Q.Cell	916144	/1pcs

Standards



SPECTRO-96, Fluorescent Spectrophotometer

The measuring sensitivity of fluorescent spectrophotometry is 2-3 order higher than that of UV VIS spectrophotometry. It can offer the information of excitation & emission spectrum, light emitting intensity & life, quantum yield rate, fluorescence polarization etc. and a wide linear range of operating curve, so that it has become an important means involved in trace analysis. It can be applied to the fields below: Medical & clinical tests—Clinical analysis of biogenic testing material • Pharmaceuticals & pharmacology—Natural pharmacon analysis, pharmaceutical

quality control and pharmacon metabolism research • Biochemistry—Minute quantity substance measurement in biological body • Food industry—Measurement of trace components in food • Pollutant analysis—Air pollution, environmental hygiene detection, food pollutants etc. • Organic & inorganic chemistry—Analysis of trace components which can not be measured with absorptive spectrophotometry.

Features: The measuring of 1200 line grating, large aperture aspherical reflector, very high sensitivity • High stable & durable xenon lamp and its power supply assure high stability of measurement and wide spectral range • Photoelectric multiplier of high performance can cover overall UV VIS region and has optimal SIN ratio • Auto null adjustment and background subtraction • Operating prompt with menu, arbitrarily setting up scale, adopting real time display of

spectrogram in sequential or compound manner (96CRT) • Spectrogram access can be implemented at any time and complex spectral processing & arithmetic operation between spectrum and data can also be implemented (96CRT) • JAnalytical conditions, spectrogram and quantitative data can be stored or read-out on hard (or floppy) disc (96CRT).

Specifications:

Model	SPECTRO-96
Wavelength	Ex200~750nm (optional optical filter) *Em200~800nm (C-T monochromator)
Wavelength Accuracy	EM±2nm
Wavelength Repeatability	EM±1nm
Spectral Slit Width	EM10nm
Sensitivity	Water Raman's peak detection with 1cm crystal fluorescent cuvette S/N≥50 (10nm narrow peak)

Features: New type cube-corner Michelson interferometer features smaller size & more compact structure, providing higher stability and less sensitive to vibrations and thermal variations than conventional Michelson interferometer • Fully sealed damp and dust proof interferometer, adopting high performance, long lifetime sealing material and desiccator, ensures higher adaptability to the environment and increases accuracy

and reliability in operation. Viewable window for silica gel enables easy observation and replacement • Isolated IR source and large space heat dissipation chamber design provides higher thermal stability. Stable interference is obtained without the need of dynamic adjustment. • High intensity IR source adopts a reflex sphere to obtain even and stable IR radiation • Cooling fan stretch suspending design ensures good mechanical stability • Super wide sample compartment provides more flexibility to accommodate various accessories • The application of programmable gain amplifier, high accuracy ND converter and embedded computer improves the accuracy and reliability of the whole system • The spectrometer connects to PC via a USB port for automatic control and data communication, fully realizing plug-and-play operation • Compatible PC control with user friendly, rich function software enables easy,

convenient and flexible operation. Spectrum collection, spectrum conversion, spectrum processing, spectrum analyzing, and spectrum output function etc. can be performed • Various special IR libraries are available for routine search. Users can also add and maintain the libraries or set up new libraries by themselves • Accessories such as Defused/Specular Reflection, ATR, Liquid cell, Gas cell, and IR microscope etc can be mounted in the sample compartment.

FTIR510 FTIR Spectrophotometer



Model	FTIR510
Spectral range	7800 to 350 cm ⁻¹
Resolution	Better than 0.85 cm ⁻¹
Wavenumber precision	±0.01 cm ⁻¹
Scanning speed	5 step adjustable for different applications
Signal to noise ratio	Better than 15,000:1 (RMS value, at 2100 cm ⁻¹ , resolution: 4 cm ⁻¹ , detector: DTGS, 1min. data collection)
Beam splitter	Ge coated KBr
Infrared source	Air cooled, high efficiency, reflex sphere module
Detector	DTGS
Data system	Compatible computer
Software	FT-IR software contains all routines needed for basic spectrometer operations, including library search, quantitation and spectrum export
IR library	11 IR libraries included
Dimensions	54x52x26cm
Weight	28kg



SPECTRO-300X, Portable FTIR Spectrometer

Model	SPECTRO-300X
Wavelength range	7000 to 400 cm^{-1} (without ATR accessory) 6000 to 600 cm^{-1} (with ZnSe ATR accessory)
Resolution, standard	1 cm^{-1}
Resolution, option	0.5 cm^{-1}
Sample compartment	Integrated
Accessories	ATR, reflection, transmission
Operating system	Windows based
Power	90-230 VAC, 12 VDC, 40 W
Dimensions (cm)	W49xD39xH20
Weight	18kg
Temp. environment	15 – 28°C
Humidity environment	Best below 65%

General:

The Interspec SPECTRO-300X series of FTIR portable spectrometers represent a low cost Fourier transform infrared and near infrared portable spectrometers and employ a number of unique features that ensure high performance from a compact instrument. The Interspec SPECTRO-300X measures just 49 X 39 X 20 cm and is regarded as one of compact and versatile infrared FTIR spectrophotometers.

The design of the SPECTRO-300X is unique both in terms of optical design and the software and firmware designed specifically to significantly reduce overall analytical times. The interferometer geometry is employing a new compact Michelson self compensating optical system that eliminates many of the optical alignment problems found in conventional type optical interferometers.

The Interspec SPECTRO-300X design avoids the use of conventional corner cube optics and dynamic alignment. In practice this means that the instrument can be used in the research laboratory, in any university or college environment and if required, can also be used outside laboratory or in remote locations.

Interferometer Performance:

All Interspec FTIR instruments offer high S:N ratios and can provide SNR up to 12000:1. Resolution in the infrared is available 2 cm^{-1} and programmable up to 32 cm^{-1} (option 0.5 and 1 cm^{-1}).

The overall wavelength range is 7000 to 400 cm^{-1} (IR) or 15000 to 3850 cm^{-1} (NIR).

The Sample Compartment:

There is no sample compartment in classical sense.

Extending Wavelength Ranges:

In order to facilitate the use of more than one beam splitter or detector, provision has been made to interchange the beam splitter and detector assemblies allowing the Interspec SPECTRO-300X to be used at any wavelength from 15000 to 400 cm^{-1} .



Beam Splitters	Range subject to coatings
KBr	7,000 to 400 cm^{-1}
ZnSe	5,000 to 500 cm^{-1}
CaF ₂	10,000 to 1,000 cm^{-1}
Quartz	15,000 to 3,000 cm^{-1}

Detector Options:

The standard IR detector is a selected high sensitivity DLATGS pyroelectric design providing the highest possible signal to noise for all but the most demanding applications. In case of NIR spectral region two types of photodiodes are available: Si and InGaAs.

Software:

Interspec for Windows software is supplied on CD and provided with each system shipped. The software includes features for all standard analytical requirements including manipulation of spectral data, instrument control, plot with preview on the screen plus many others. Also included are several facilities for analytical modelling of interferograms or spectra, with smoothing, and baseline correction, interactive editing and data manipulation. Also spectral subtraction, mixture subtraction, smoothing derivatives, plot with preview etc. Data input and output is possible in ASCII or JCAMP. Other commercial programmes can be used including Thermo/Galactic GRAMS for features such as Library Search. The Interspec for Windows programme is written in 32 bit protected mode. Our unique software has been designed specifically for multi function applications, it is easy to use and it is provided free of charge. The utility of the Interspec for Windows programme can be extended by adding other commercial programmes such as search, component identification, Kramers Kronig Transform, Chemometrics, etc. to suit individual requirements.



AAS-210, Atomic Absorption Spectrophotometer

**Features:****Innovated rich oxygen air-acetylene flame analysis technique:**

The patented flame analysis technique adopting rich oxygen air-acetylene flame as the substitution for nitrous oxide-acetylene flame for high temperature element analyses, such as Ca, Al, Ba, W, Mo, Ti, V, etc. Flame temperature is continuously adjustable between 2300-2950°C, which makes it possible to choose the best atomization temperature for different elements. It features easy operation, low analysis cost and wide flame AAS analytical range. Rich oxygen flame will not pollute the environment and is not harmful to human bodies. It's a break-through in flame AAS analysis.

Integrated flame/graphite furnace atomization system, changeable with flame emission burner:

- Automatically controlled changeover of the integrated flame and graphite furnace atomizer featuring easy operation and time saving eliminates human labor.
- A flame emission burner head can be installed to perform flame emission analysis to alkali metals as K, Na etc.

Accurate fully automated control system:

- Automatic 6-lamp turret, automatic adjustment of lamp current and optimization of light beam position.
- Automatic wavelength scanning and peak picking.
- Automatic spectral bandwidth changing.
- Automatic changeover between flame and graphite furnace operation, automatic optimization of position parameters, automatic ignition and automatic gas flow setting.

Reliable fully automatic graphite furnace analysis:

- Adopting FUZZY-PID and dual curve mode light-controlled temperature control technique, temperature auto-correction technique, ensures fast heating, good temp. reproducibility & high analytical sensitivity. The temperature control accuracy is less than 1%.

- Graphite furnace with pneumatic control and pressure lock ensures constant pressure and reliable contact.
- Multi-function autosampler features automatic standard sample preparation, automatic correction of sampling probe depth, automatic tracing and correction of liquid surface height in the sample vessel, with the sampling accuracy of 1% and reproducibility of 0.3%, realizing fully automation of graphite furnace analysis.

Perfect safety protection measures:

- Alarm and automatic protection to fuel gas leakage, abnormal flow, insufficient air pressure and abnormal flame extinction in flame system.
- Alarm and protection function to insufficient carrier gas and protective gas pressure, insufficient cooling water supply and over-heating in graphite furnace system.

Advanced and reliable electronic design:

- Adopting large-scale programmable logic array and Inter I2C bus technology.
- European type sockets and AMP adapters with high reliability to ensure long term reliability of the whole electronic system.

Easy and practical analysis software:

- Easy-to-use MS analysis software is made under Windows operating system, realizing fast parameter setting and optimization.
- Automatic sample dilution, automatic curve fitting, automatic sensitivity correction.
- Automatic calculation of sample concentration (content), mean value, standard deviation and relative standard deviation calculation.
- Multi-elements determination in sequence to the same sample.
- Measured data and final results can be printed out and edited in Excel format.

For detail information look at www.mrclab.com